

## Certificate in relation to air transportation by passengers with walking disabilities

We hereby certify that our products

### **e-motion M10 and e-motion M12**

including the nickel metal hydride batteries (Special Provision A123) fulfil the requirements of

### **IATA chapter 2.3.2.2**

and therefore may be carried together with the wheelchair as **checked baggage** without restriction.

Condition precedent to meet the above listed requirement is the transportation in a way that an unintentionally switching on and damage by other luggage is avoided. The batteries can be removed for transportation (protected from damage and from short circuit) and must be carried in in strong, rigid packaging in the cargo compartment.

The pilot-in-command must be informed of the location of the battery.

**We strongly recommend obtaining the approval by corresponding Airline Company before commencement date respectively departure date.**

This certificate does not apply to damaged or defect batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31<sup>st</sup> December 2025**.

Albstadt, 20<sup>th</sup> of November 2024

Alber GmbH



Johannes Beltman  
Managing Director