



Certificate in relation to air transportation by passengers with walking disabilities

We hereby certify that our product

smoov (one) O10

including the built-in lithium-ion batteries fulfil the requirements of

IATA chapter 2.3.2.4 (Wheelchair / Mobility Aids)

The lithium-ion batteries are tested according to UN Manual of Tests and Criteria, Part III, Subsection 38.3.

For lithium batteries installed in mobility aids, there is no Watt-hour limit.

Additional information: The rated energy of the built-in drive unit battery (not removable by the user) is **below 300 Wh**. The rated energy of the control unit (item number 1594326) battery (also not removable by the user) is **below 10 Wh**.

Since drive unit and control unit are portable and contain batteries, we recommend to ask, if you can take them with you to the passenger cabin. Both units must be switched off before boarding and protected from damage and unintentional activation.

The pilot in command must be informed about the location of the mobility aid and the control unit.

We strongly recommend contacting the airline sufficiently in advance to request approval for transport.

This certificate does not apply to damaged or defective devices or batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31st December 2025**.

Albstadt, 22nd of November 2024

Alber GmbH

A handwritten signature in blue ink, appearing to read "J. Beltman", with a long horizontal flourish extending to the right.

Johannes Beltman
Managing Director